


ULTIMA WEB REMOTE ACCESS

Share secure data and images


Ultima's Web Remote Access (WRA) interface provides access to Ultima Terabox studies from outside the facility. Because protecting these data is essential, the solution offers advanced user account management as well as secure communication.



The login interface is titled "Ultima Web Remote Access". It features a white background with a dark blue border. There are two input fields: "Username:" and "Password:". Below the password field is a "Login" button.

Set user's rights according to specific rules:

- Authorized users (patients, doctors, etc.)
- Temporary or permanent authorization
- User's accounts rights (viewing images, reports, old studies, etc.)
- Sort by types of studies (scanner, X-ray, MRI, etc.)




The "User Configuration" window is divided into two main sections. The left section, titled "Account Information", contains fields for "User Name:" (with value "user"), "Password:", "Password2:", "Group:" (with a dropdown menu showing "Patients (Patients)"), and "Type:" (with radio buttons for "Administrator" and "User", where "User" is selected). The right section, titled "Linked to:", contains radio buttons for "Patient" (selected), "Performing Physician", "Referring Physician", and "None". Below these is a "Linked To values:" field with a dropdown menu and "Add" and "Remove" buttons. A list box below shows the value "pat0010". At the bottom of the window are "Save" and "Cancel" buttons.

Key Features

- Searches studies in the system
- Displays the study report
- Displays old patient studies
- Integrated web viewer
- Transmits secure data
- Configurable compression rate before sending
- Limited access to view studies
- Limited access to view studies' reports
- Limited access to create and edit studies' reports
- Interface assignment and dismissal of temporary rights for system administrators
- Advanced user's rights
- Complete audit log of actions performed by system users

Web Viewer

The Ultima Web Viewer software is integrated into the system. Its user-friendly, fast and ergonomic interface includes a full range of measures and annotations, as well as an advanced comparative mode.



MED SQUARE Ultima Web Remote Access									
Users Groups Permissions Connections Log admin Logout									
Patient ID	Patient Name	DOB	Acc. Number	Description	Modality	From	To	Compression	
					All			JPEG Lossless	
								Search	
Patient ID	Patient Name	Description	Date/Time	Modalities	Description	Date/Time	Modalities		
14221420110211	TEST		2/11/2011 2:22 PM	USIR					
432863	20101102110211		2/11/2011 11:41 AM	USIR					
127643	20101102110211		2/11/2011 11:12 AM	USIR					
279006	20101102110211		2/11/2011 10:28 AM	USIR					
411909	20101102110211		2/11/2011 10:05 AM	USIR					
432852	20101102110211		2/11/2011 9:26 AM	USIR					
389082	20101102110211		2/11/2011 8:55 AM	USIR					
427913	20101102110211		2/10/2011 6:16 PM	USIR					
12	20101102110211	TSA CR RL	2/9/2011 6:14 PM	CT					
12	20101102110211	TSA CR RL	2/9/2011 6:11 PM	CT					
20110209115002	20110209115002	TOTO	2/9/2011 11:51 AM	CT					
PKR730000	20110209115002		1/20/2011 3:16 PM	NM					
ALGO00001	20110209115002		1/20/2011 3:16 PM	CT					
FUJ000006	20110209115002		1/20/2011 3:01 PM	CR					
FUJ000001	20110209115002		1/20/2011 3:01 PM	CR					

Medsquare, 17 rue du Jura, 75013 Paris, FRANCE

W medsquare.com E contact@medsquare.com T + 33 (0)1 55 25 62 50