


# ULTIMA WEB VIEWER & ULTIMA CD VIEWER

## At the centre of your images distribution

### Ultima Web Viewer **Ultima CD Viewer**

These two viewers have the same interface and functionalities, but they differ in their medium of distribution:

- **Ultima Web Viewer** enables you to distribute a complete study (images and reports) inside (intranet) and outside (internet) your facility through **web pages**.
- **Ultima CD Viewer** enables you to distribute a complete study (images and reports) inside (intranet) and outside (internet) your facility through **CD, DVD or Blu-Ray**.



#### Common key features:

- Multi-modalities viewer: Single frame (CT, MR, CR, etc.) and multi-frame (US, XA, etc.)
- Multilingual interface: English, French, Spanish, Italian, German, Portuguese
- Viewer starts automatically without any installation (work on a secured PC)
- View of the report of the on-going study
- Key object selection and presentation state



### So versatile and easy to use:

- Supports images from any modality (CT, MR, XA, US, etc.)
- Multi-monitor support
- Windowing mode pre-set by type of modality and accessible at anytime
- Hanging protocols module for optimal diagnostics
- Viewing window fully customizable (number of series per screen, number of images per series, etc.)
- Full range of measurements and annotations (distance, angle, text, ROI, etc.)
- Simplified creation of boards of images for printing on a Windows printer
- Evolving comparative mode (intra and inter studies) in automatic (based on the positioning of the image in space) or manual (ability to create groups of series)
- Viewing and advanced processing feature: CT/MR automatic reference line and XA (subtraction angiography)
- Multi-planar reconstruction (MPR)



Medsquare, 17 rue du Jura, 75013 Paris, FRANCE

**W** [medsquare.com](http://medsquare.com) **E** [contact@medsquare.com](mailto:contact@medsquare.com) **T** + 33 (0)1 55 25 62 50

© 2017 Medsquare SAS — All rights reserved.

Medsquare SAS reserves the right to modify the design and specifications described in this document at any time without notice.

FP17us/01-jun17